

Instructions for Use H2O Cem



IMPORTANT

Read this 'Instructions for Use' document before handling. This product is for use by a dental professional only.

Adheres chemically to tooth substance and non-precious alloys to optimise retention. Rapidly increasing pH during setting produces excellent biocompatibility. Contains fluoride ions. Superb translucency for optimum aesthetics under glass-ceramic and porcelain jacket crowns.

IDENTIFICATION

Product Name	H2O Cem
Company Name	Perfection Plus Ltd
Company Address	6 Westwood Court, Totton, Hampshire, SO40 3WX, UK
Company Phone No.	+44 (0) 2380 866 677

PRODUCT CODE

REFPP010/1550 H2O Cem 30g Powder

INDICATIONS FOR USE

For cementing crowns, bridges, inlays, onlays and orthodontic bands.

INSTRUCTIONS FOR USE

Mixing: The recommended powder: liquid ratio is 3.6:1 mm at 22-24°C and RH 40-60%. Use a clean and dry polished glass mixing slab and a stainless steel flat-bladed spatula. Invert bottle to "fluff" powder for accurate dispensing. Measure ONE scoop of powder onto glass slab, taking care not to compress powder against side of bottle with the scoop. Remove excess powder from the scoop using a flat-bladed spatula. Dispense TWO drops of water onto slab. To deliver accurate water drops, hold bottle vertically. If water does not drop from the "tip" of steel tube, clean with a tissue. Divide powder into two equal portions and incorporate half the powder into water and mix for 10-15 seconds, then add remaining powder and spatulate to a creamy consistency. Total mixing time is 30 seconds. DO NOT ADD POWDER IN SMALL INCREMENTS.

Working Time: A minimum of 1 minute 45 seconds from start of mix at 22-24°C.

Clinical Setting Time: A maximum of 6 minutes from end of mix.

Placement: Apply the cement on completion of mix, using normal instruments. Excess cement may be removed immediately or after cement has hardened. DO NOT ATTEMPT TO REMOVE CEMENT WHILST IN RUBBERY PHASE AS ADHESION MAY BE IMPAIRED. Clean instruments immediately on completion.

WARNING!

Avoid contact of powder and mixed cement with eyes. In case of contact, wash thoroughly with water and obtain medical advice.

STORAGE

Store powder in a cool, dry place. The powder is protected against moisture uptake by a desiccant capsule, since moisture ingress can adversely affect the shelf life. Always replace caps immediately after use.



Batch code/Expiry Date

The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.

The expiry date is shown in year, month format. Do not use the product after this date.

TESTING

For test purposes, the powder: liquid ratio is: 1.08g to 0.30g at 23±1°C RH 50±10%. ISO 9917 net setting time: A maximum of 5 minutes from end of mix at 36-38°C.

SPECIALLY FORMULATED FOR USE IN DENTISTRY

