

Latex Dam Instructions for Use

IMPORTANT

Read this 'Instructions for Use' document before handling. This product is for use by a professional only.

SPECIALLY FORMULATED FOR USE IN DENTISTRY

IDENTIFICATION

Product Name	Latex Dam 6 x 6 pk 36
Company Name	Perfection Plus Ltd
Company Address	6 Westwood Court, Brunel Road, Totton, Hampshire, SO40 3WX, UK
Company Phone No.	+44 (0) 2380 866 677

PRODUCT CODE – REF

004/0032 Latex Dam 6 x 6 Thin Green pk 36

004/0033 Latex Dam 6 x 6 Medium Green pk 36

004/0034 Latex Dam 6 x 6 Heavy Green pk 36

INDICATIONS FOR USE

Latex rubber dams are manufactured from natural rubber latex. The dental dam is used as a protective oral barrier during endodontic procedures to reduce the risk of infection, to isolate the working area, to prevent particles or small objects making their way into the patient's airways and to improve patient comfort. The dental dam is a single use, non-invasive, non-sterile device that is designed to be used only by qualified dental professionals for the described purpose.

CONTRAINDICATIONS

It is the responsibility of the dental practitioner to determine if the use of dental dam is appropriate for each individual patient. Latex dental dams should not be used in situations where the skin in the relevant area is broken, the patient has a known latex allergy or sensitivity to type IV allergens or other clinical considerations indicate that it should not be used.

INSTRUCTIONS FOR USE

1. Take a sheet of latex dam and a sheet of rubber dam template.
2. Place latex dam on top of template and adjust.
3. Using a marker pen, select which tooth is required for isolation.
4. Make a hole using Perfection Plus's dental dam punch.
5. Stretch latex dam on to a metal dam frame.
6. Place dental dam clamp in the hole and position around the tooth.

The dental dam punch must be in good condition, as failure to punch a clean hole could result in the dental dam tearing.

STORAGE


The product is flammable and should be kept away from naked flames during handling and storage.

Store below 26°C and away from direct sources of heat.

Avoid damp storage conditions and keep out of direct sunlight.

Avoid contact with copper, copper containing alloys, oils, solvents and greases.

The box should be closed when not in use.



WARNING!

Natural rubber latex dental dams should not be used on persons who are susceptible to natural rubber latex allergies but is otherwise considered non-hazardous and safe for use.

The product contains processing additives including dithiocarbonate that may cause sensitisation by skin contact and give rise to Type IV Allergies. Persons who know they are susceptible to dithiocarbonates should avoid contact with the product.

Talc is considered non-hazardous in the quantities present on the product, but contact should be avoided with broken skin.

The product is only to be used by a qualified dental practitioner trained in the use of dental dams.



- Do not use the dental dam pad if there is any sign of damage. Please dispose of any damaged items according to local regulations
- Do not use the dental dam product for any other purpose that that stated in this manual.
- The product should only be used before the expiry date indicated on the box.
- Do not use the product on children unless under the supervision of a qualified practitioner.
- This product is for single use only. It is not suitable for cleaning

Caution



Should a serious incident in relation to the product occur, it should be reported to the manufacturer and the competent authority.

LOT NUMBER AND EXPIRY DATE

The above are printed on each package and should be quoted in all correspondence which requires identification of the product. Do not use the product after this date. See product box for lot.



Perfection Plus Ltd
6 Westwood Court
Brunel Road
Hampshire
SO40 3WX, UK
www.perfectionplus.com



Perfection Plus EU Ltd
The Black Church
St Mary's Place.
Dublin, Ireland,
D07 P4AX
Regulatory@perfectionplus.com



Advena Ltd
Tower Business
Centre
2nd Floor Tower
Street
Swatar, Malta

