



## PoliFLEX Plus

### Instructions for Use

#### **IMPORTANT**


Read this 'Instructions for Use' document before handling. This product is for use by dental professionals only.

#### **SPECIALLY FORMULATED FOR USE IN DENTISTRY**

#### **IDENTIFICATION**

The Instructions for Use are for single use medical devices and only covers the products listed in the relevant product group schedules.

The Class 2a medical devices are for use in the mouth by or under the instruction of a qualified dental professional and are manufactured in accordance with the standard ISO 1797, and should only be used in conjunction with a rotary hand piece that conforms to ISO 14457. The operator should ensure that the rotary instrument is correctly installed in the hand piece prior to commencement of any procedure.

Product Name	PoliFlex Plus
Company Name	Perfection Plus Ltd
Company Address	 6 Westwood Court, Brunel Road, Totton, Hampshire, SO40 3WX, UK
Company Phone No.	+44 (0) 2380 866 677

#### **PRODUCT CODE – REF**

005/5000	Kit 240 Assorted Discs & 2 Mandrels	005/5007	Refill 14mm Coarse Pk 100
005/5001	Refill 10mm Super Coarse Pk 100	005/5008	Refill 14mm Medium Pk 100
005/5002	Refill 10mm Coarse Pk 100	005/5009	Refill 14mm Fine Pk 100
005/5003	Refill 10mm Medium Pk 100	005/5010	Refill 14mm Ultra-Fine Pk 100
005/5004	Refill 10mm Fine Pk 100	005/5011	Mandrels Pk 10 RA
005/5005	Refill 10mm Ultra- Fine Pk 100	005/5012	Mandrels Single RA
005/5006	Refill 14mm Super Coarse Pk 100		

#### **INDICATIONS FOR USE**

PoliFLEX Plus discs when attached to RA mandrels are for used for finishing and polishing of composite restorations, gold and amalgam. The discs are produced in 2 different shapes and 5 different grits and give best performance when applied at slow speed (either at maximum of 10,000 or 30,000 RPM dependant on the grit size), under slight pressure.

In accordance with the Medical Device Directive 93/42/EEC and Medical Device Regulation 2017/745, these devices are invasive through the body orifice not beyond the larynx, are not implantable and are intended for transient use. The maximum number of repeat applications required is determined by the qualified dental professional. The area of contact within the mouth is the teeth.

There are no claims made by the manufacturer other than for the cleaning of teeth and improved dental hygiene.

#### **WARNINGS/CONTRAINDICATIONS**

Used rotary instruments should be considered as contaminated and appropriate handling precautions should be taken. Gloves, eye protection and a mask should be worn. Other measures may be required if there are specific infection or cross-contamination risks from the patient.

These devices contain small amounts of nickel and should, therefore, not be used on individuals with a known sensitivity to this metal as in extreme cases it may cause hypersensitivity.

**STORAGE (pre-use):** The instruments should be stored in the container (stand or pouch) until required. Containers or pouches must be dry before opening to avoid recontamination of the contents with water. Storage should be in dry, clean conditions and at ambient temperature.

#### **STERILISATION**

The discs do not require sterilisation. The rotary instrument should be kept in its original packaging at room temperature and protected against dust and moisture until used for the first time.

#### **INSPECTION**

Inspect the instruments, with the aid of magnification if necessary, and discard any damaged or corroded instruments.

## **PROPER USE**

- The instruments are for single use only. Do not reuse.
- Inspect the device prior to use for any defects.
- Only use hand piece, angles and turbines that are technically and hygienically flawless, maintained and cleaned.
- Ensure handpieces are maintained in good working order and remain adequately lubricated at all times to ensure maximum effectiveness of the device. Failure to properly maintain the handpiece can lead to procedural delays, injury to the user or patient, aspiration or swallowing of the device or damage to the preparation site.
- The device and hand piece must be concentric and true running. Instruments that are deformed or no longer run true should not be used and must be disposed of.
- The instruments must be fully inserted into the handpiece and locked where applicable.
- The instruments are to be brought to speed before placing on the object.
- Polish with gentle circular movements.
- Avoid tilting or levering because of the increased risk of breakage.
- Any damaged and incorrectly shaped instruments will cause vibration.
- Instruments that are deformed or no longer run true should not be used and must be disposed of.
- Wear eye protection.
- Wear a respiratory mask to prevent inhaling any dust and/or debris during the procedure.
- These products must only be used by qualified dental professional.
- Incorrect use produces poor results and increases risk.
- Incorrect use may harm tissue, cause premature wear, destroy the instrument and endanger the operator, patient or third parties.
- The devices are designed and manufactured to perform safely when used in combination with CE marked medical devices (i.e. prophylaxis paste or mandrels if required) which are themselves intended for use in the oral cavity.
- Do not use the discs at speeds in excess of the recommended RPM.
- Do not operate the mandrel without the disc attached
- Avoid touching composite with the mandrel or disc eyelet because discolouration may occur. This discolouration can be removed by repetition of the polishing steps. The smaller, low profile, mandrel head and unique eyelet design reduces this risk of composite and eyelet contact.
- Use the discs in the correct order, incorrect procedures can result in a reduction in the polishing quality.
- To ensure vibration-free working the connected instrument must be centred after mounting on the mandrel or shaft. Any damaged and incorrectly shaped instruments will cause vibration. Bent or non-concentric instruments must be disposed of.
- PoliFLEX Plus discs are to be used with the PoliFLEX Plus Mandrel or equivalent. PoliFLEX Plus mandrels and polishing discs are proven to function together correctly by historical use with no adverse reports.

## **VALIDATION**

These processes have been validated as being capable of preparing Stoddard dental instruments for use. It remains the responsibility of the processor to ensure that the processing as performed using the equipment, materials and personnel achieve the required results. This may require validation and monitoring of the process. Any deviation from these instructions should be properly evaluated for effectiveness and potential adverse results.

## **PRESSURE**

Avoid excessive pressing force!

- Excessive pressure must be avoided at all times.
- Excessive pressure may damage the working sections of rotary polishing instruments or damage the fill material. Heat build-up is increased

## **INSTRUCTIONS/RECOMMENED SPEEDS**

<b>Grit Size</b>	<b>Maximum RPM</b>
Super Coarse	10,000 rpm
Coarse	10,000 rpm
Medium	10,000 rpm
Fine	30,000 rpm
Ultra-fine	30,000 rpm

- 1) Always use PoliFLEX Plus discs when attached to a mandrel with a slow speed hand piece (maximum 10,000 or 30,000 RPM). Affix the disc to the mandrel by gently pushing the eyelet piece on to the mandrel until the disc is securely in place.
- 2) The polishing motion should be constant and uni-directional, e.g. commencing at the gingiva and moving outwards over the restoration. A back and forth movement over the composite – enamel margin is not recommended.

- 3) Use light pressure when polishing, let the discs do the work.
- 4) Keep the restoration surface and disc dry while polishing. A dry surface will produce a smoother, more uniform surface.
- 5) Remove PoliFLEX Plus discs from the mandrel by:
  - a) positioning a thumbnail under the disc eyelet portion and pushing the disc away from the hand piece (i.e. Pop off) or;
  - b) grasping the disc and eyelet and peeling the disc upwards away from the handpiece.

#### **INTENDED PATIENT GROUPS/INTENDED PURPOSE**

The devices are intended for use on patients of any age. The products are only aimed as a tool for use by the qualified dental professional who is responsible for determining the treatment required by each individual patient and which treatment would outweigh the risks of performing.

The final brush selection will be dependent on 3 factors:

- the dental professional performing the procedure,
- the patient undergoing treatment,
- the type of procedure being performed.

#### **INTENDED USERS**

All devices are intended to be used by a qualified dental professional.

#### **CAUTION**

Should a serious incident in relation to the product occur, it should be reported to the manufacturer and the competent authority.



#### **BATCH CODE/EXPIRY**

The batch code is shown as LOT. Please quote this in all correspondence.

The product lifetime is estimated to be 3 years from the date of manufacture but may vary depending on the method of use (cutting, conditions), the frequency of use, or when the sterilisation has reached a maximum of 10 cycles.

#### **LIFETIME OF THE DEVICE**

The lifetime of the device is normally determined by wear and tear and depends on a number of factors and actions being carried out by the end user including:-


- the instructions for use being followed correctly to avoid damage occurring to the instrument and/or patient,
- the condition of the handpiece being used,
- possible solutions being used with the product,
- the device being routinely inspected for defects before any procedure

Stoddard products are tested and certified as both safe and usable within the date parameters set by the use by dates as indicated on the product labelling. Whilst we would not expect any real degradation of the product, we cannot guarantee its safety beyond the use by date and recommend that the products are disposed of in accordance with the local regulations. The use by date set for these devices is 5 years from the date of manufacture. The 5 year shelf life has been concluded from the inspection tests for dimensional checks, shank/mandrel critical connection checks and shelf life tests conducted on products over 5 years old and the product performing correctly and without failure demonstrating the 5 year shelf life is acceptable.


Stoddard products are manufactured and packaged in such a way that no deterioration can occur which could compromise the safety of the patient, user or other persons or the performance of the product.

#### **SAFETY PRECAUTIONS**

These dental instruments were developed and manufactured for their specific dental surgical application. Incorrect use may harm tissue, cause premature wear, destroy the instrument and endanger the operator, patient or third parties.

 Stoddard  
Manufacturing  
Company Ltd  
Black Horse Road  
Letchworth  
Herts, SG6 1HB, UK

 Perfection Plus Ltd  
6 Westwood Court  
Brunel Road, Totton  
Hampshire  
SO40 3WX, UK  
[www.perfectionplus.com](http://www.perfectionplus.com)

 EC REP  
Stoddard ApS  
Moellegade 32  
8000 Aarhus C  
Denmark

 Perfection Plus EU Ltd  
The Black Church  
St Mary's Place.  
Dublin, Ireland,  
D07 P4AX  
Regulatory@perfectionplus.com

 MD  
 CE  
1639