

Instructions For Use PerfectEtch 12g



IMPORTANT

Read this 'Instructions for Use' document before handling. This product is for professional use only.

IDENTIFICATION

| | |
|-------------------|----------------------------------------------------------------|
| Product Name | PerfectEtch 12g and 20 tips |
| Company Name | Perfection Plus Ltd |
| Company Address | 6 Westwood Court, Brunel Road, Totton, Hampshire, SO40 3WX, UK |
| Company Phone No. | +44 (0) 2380 866 677 |

PRODUCT CODE

| |
|-----------------------|
| REF PP010/1001 |
|-----------------------|

INDICATIONS FOR USE

| |
|-------------------------------------|
| Acid etching of enamel and dentine. |
|-------------------------------------|

CONTRAINDICATIONS

| |
|----------------------------------------------------------------------------------------------------------------|
| PerfectEtch 12g are contraindicated for use on soft tissue including mucosa, skin, eye and dental pulp tissue. |
|----------------------------------------------------------------------------------------------------------------|

INSTRUCTIONS FOR USE

| |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ol style="list-style-type: none"> 1. Remove the cap from the syringe and save. Attach a new applicator tip provided. 2. Prior to intraoral application, dispense a small amount of etch onto a mixing pad to verify flow and ensure the applicator tip is properly secured. 3. Prepare tooth with plain pumice or oil-free and fluoride-free prophylaxis paste, such as Perfection Plus Oil Free Prophylaxis Paste, in the usual manner. 4. Rinse thoroughly, isolate and dry the tooth with oil-free compressed air. 5. Apply etchant to the tooth surface for 15 seconds. If you are applying to dentine, etch for 10 seconds. Never leave the etchant gel on for more than 60 seconds. Lightly agitating the etchant during its application can enhance its action. 6. Rinse thoroughly with water. The dry etched area should appear dull and frosty white. If this is not the case, repeat step 6. |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

WARNING!

| |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>WARNING! PerfectEtch 12g contains phosphoric acid which is corrosive. Avoid contact with soft tissue, eyes and skin. If contact occurs flush immediately with water for at least 15 minutes. Seek immediate medical attention.</p> |
| |

PRECAUTIONS

| |
|----------------------------------------------------------------------------------------------------------------------------|
| Wear proper personal protection equipment before using PerfectEtch 12g, to prevent contact with soft tissue, eyes or skin. |
|----------------------------------------------------------------------------------------------------------------------------|

ADVERSE REACTIONS

| |
|-----------------------------------------------------------------------------------------------------------------|
| Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to PerfectEtch 12g. |
|-----------------------------------------------------------------------------------------------------------------|

STORAGE

| |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Recommended storage conditions are between 2°C to 27°C, bringing to room temperature before using. (23°C). Do not freeze. Always replace cap after each use.</p> |
| |

Batch code/Expiry Date

| |
|-------------------------------------------------------------------------------------------------------------------------------|
| The batch number should be quoted in all correspondence. See packaging for lot and expiry date. Do not use after expiry date. |
|-------------------------------------------------------------------------------------------------------------------------------|

SPECIALLY FORMULATED FOR USE IN DENTISTRY

