

## Instructions for Use Saliva Ejectors



### **IMPORTANT**

Read this 'Instructions for Use' document before handling. This product is for professional use only.

### **IDENTIFICATION**

Product Name	Saliva Ejectors
Company Name	Perfection Plus Ltd
Company Address	 6, Westwood Court, Totton, Hants. SO40 3WX UK
Company Phone No.	+44 (0) 2380 866 677

### **PRODUCT CODE**

**REF** PP011/4250 Protect+ Saliva Ejectors Clear 150mm Bag 100

**REF** PP004/4251 Protect+ Saliva Ejectors White 150mm Bag 100

**REF** PP004/4252 Protect+ Saliva Ejectors Blue 150mm Bag 100

### **INDICATIONS FOR USE**

For use by dental professionals to eject saliva from a patient's mouth.

### **INSTRUCTIONS FOR USE**

Insert firmly into the conic connector of aspiration tube.  
Single Use Only.

### **WARNING**

Do not use if tip is loose.  
Professional use only.



### **Batch code/Expiry date**

The batch number should be quoted in all correspondence. See packaging. The useful life of the product is 5 years but there is no expiry date.

**SPECIALLY FORMULATED FOR USE IN DENTISTRY**

