

ImpressPLUS-5

5 day Extra Dimensional Stability Alginate



ImpressPLUS-5 is a non-chromatic alginate with extra dimensional stability of 120 hours.

The orange colour allows an excellent contrast and aids the ability to evaluate the quality of the Impression.

The impression can be disinfected by submersion in specific solutions. Follow the manufacturer's instruction for use to ensure no alteration.

Mixing Times	
Mixing Time with spatula	45 seconds
Working Time to load the impression tray	1 minute 45 seconds (including 45 seconds Mixing Time)
Initial Setting Time for inserting in the mouth	1 minute 55 seconds (including 1 minute 45 seconds Working Time). Only 30 seconds in the mouth.

Delivery System

- 1x 450g bag

Presentation

REF 010/1299

Characteristics

- Class A type 1 and type 2 alginate
- Easy to mix
- Thixotropic
- Mint flavour
- Dust Free
- 5 years expiry date

Technical Data

- Extra fast in mouth, minimum setting time 30"
- High accuracy: 20 micron
- Extra dimensional stability (120 hours)