ImpressPLUS

Alginate





ImpressPLUS-5 Alginate

A non-chromatic alginate with extra dimensional stability for 5 days. The orange colour allows excellent contrast to evaluate the quality of the impression.

Delivery System

• Bag 450g

Presentation REF 010/1299

Characteristics

- Class A
- Easy to mix
- Thixotropic
- Extra fast in mouth, minimum setting time 30 seconds
- Mint flavour
- Dust Free

Technical Data

Very high accuracy: 20 micron

Extra dimensional stability (120 hours/ 5 days)

High elasticity



ImpressPlus Chrome Alginate

Fast setting Chromatic alginate with colour changing from Violet through to Lilac then Green indicating the individual phases to guarantee the mixing and insertion times.

The short setting time (60 seconds) minimises any patient discomfort whilst the high elasticity and resistance to tear provide a high-definition impression.

ImpressPLUS

Alginate





ImpressPLUS Alginate

Mint flavoured, dust free alginate available in two settings. Fast setting (60 seconds) minimises any patient discomfort whilst the regular setting allows a longer working time to mix and load the impression tray. Both offer high elasticity and tear resistance along with high detail definition.

Delivery System

• Bag 450g

PresentationREF010/1295Fast SettingREF010/1294Regular Setting

Characteristics

- Class A
- Easy to mix
- Thixotropic
- Color after setting: green
- Mint flavour
- Dust Free

Technical Data

High accuracy: 50 micron

Alginate Accessories



Water and Alginate Measures

Ref 010/1297

Fix II Tray Adhesive







Alginate Spatulas

Plastic Spatula White
Assorted Spatulas Pack 6
Plastic Spatula Lime
Plastic Spatula Light Blue
Plastic Spatula Glitter Blue
Plastic Spatula Green
Plastic Spatula Pink
Wooden Handled Spatula Straight
Wooden Handled Spatula Curved

Perfection Plus Limited 6 Westwood Court, Brunel Road, Totton, Hampshire SO40 3WX, United Kingdom Email sales@perfectionplus.com www.perfectionplus.com